

Restriction sites for pGEX-6P-1

The following restriction site table was compiled using DNASIS software for sequence analysis. The enzymes chosen are those which we believe to have been commercially available in June 1992.

The vector pGEX-6P-1 was not tested with each enzyme and therefore the accuracy of the table cannot be guaranteed. Please contact your local Pharmacia subsidiary if a discrepancy is identified.

Please note that enzymes with non-palindromic recognition sequences are listed twice. You must combine both entries to obtain the total number of sites for these enzymes.

Locations represent the 5' end of the recognition sequence for each restriction enzyme.

Enzyme (# sites): Locations

Aat II (1): 1256
Acc I (1): 964
Acc65 I (0)
Aci I (19): 978, 1335, 1664, 1831, 2475, 2619, 3115, 3136, 3143, 3242, 3395, 3482, 3711, 4182, 4238, 4339, 4349, 4666, 4871
Aci I (39): 225, 549, 974, 1080, 1119, 1129, 1171, 1288, 1434, 1543, 1620, 1785, 2022, 2113, 2484, 2729, 2850, 2869, 2996, 3024, 3186, 3203, 3229, 3283, 3434, 3450, 3459, 3501, 3573, 3637, 3657, 3660, 3888, 3972, 4128, 4266, 4405, 4508, 4742
Afl II (0)
Afl III (3): 496, 3067, 3683
Age I (0)
Aha II (5): 34, 1256, 1638, 3639, 4322
Alu I (25): 54, 178, 482, 631, 809, 1048, 1067, 1825, 1888, 1988, 2509, 2766, 2812, 2902, 3128, 3463, 3743, 4167, 4282, 4375, 4439, 4534, 4643, 4754, 4934
Alw I (7): 914, 945, 1532, 1853, 2317, 2415, 2501
Alw I (4): 941, 1545, 2325, 2422
AlwN I (1): 2653
Apa I (2): 936, 3890
ApaL I (4): 18, 1507, 2753, 3663
Apo I (3): 954, 3248, 3958
Asc I (0)
Ase I (4): 189, 2003, 4367, 4426
Asu I (15): 290, 752, 936, 937, 1199, 1815, 2037, 2054, 2133, 3545, 3890, 3891, 4235, 4259, 4663
Asu II (1): 654
Ava I (2): 959, 969
Ava II (4): 752, 1815, 2037, 4235
Ava III (0)
Avr II (0)

Bal I (1): 463
BamH I (1): 945
Ban I (6): 2226, 3603, 4322, 4452, 4732, 4744
Ban II (2): 936, 3890
BbrP I (0)
Bbs I (3): 1189, 3822, 4161
Bbs I (0)
Bbv I (11): 42, 1046, 1143, 1935, 2124, 2726, 3145, 3163, 4373, 4622, 4695
Bbv I (7): 1746, 2440, 2646, 2649, 3732, 4103, 4232
Bcg I (1): 3986
Bcg I (2): 299, 1652
Bcl I (2): 692, 3697
Bfa I (4): 270, 1986, 2321, 2574
Bfr I (0)
Bgl I (2): 2055, 4698
Bgl II (0)
Bpm I (3): 2089, 3989, 4760
Bpm I (1): 3536
Bpu1102 I (0)
Bsa I (1): 2107
Bsa I (0)
BsaA I (1): 1159
BsaB I (0)
BsaH I (5): 34, 1256, 1638, 3639, 4322
BsaJ I (10): 932, 940, 941, 949, 959, 2907, 4318, 4457, 4594, 4884
BseA I (0)
Bsg I (2): 69, 523
Bsg I (3): 11, 3549, 3749
BsiE I (6): 975, 1658, 1807, 2730, 3154, 4672
BsiW I (0)
BsiY I (11): 264, 861, 1052, 2583, 2862, 3028, 3046, 3259, 4105, 4349, 4800
Bsl I (11): 264, 861, 1052, 2583, 2862, 3028, 3046, 3259, 4105, 4349, 4800
Bsm I (0)
Bsm I (0)
BsmA I (5): 2108, 3375, 3780, 3906, 4293
BsmA I (2): 1051, 1338
Bsp1286 I (7): 18, 936, 1507, 1592, 2753, 3663, 3890
BspD I (0)
BspE I (0)
BspH I (3): 1234, 1339, 2347
BspM I (0)
BspM I (1): 63
BspM II (0)
BspW I (26): 102, 866, 1935, 2055, 2443, 3015, 3129, 3194, 3360, 3450, 3546, 3648, 3735, 3872, 3879, 4080, 4087, 4089, 4231, 4325, 4409, 4453, 4668, 4698, 4707, 4733
Bsr I (14): 1529, 1968, 2011, 2129, 2535, 3497, 3599, 3836, 3988, 4298, 4397, 4558, 4584, 4806

Bsr I (6): 1153, 1703, 2651, 2664, 3338, 3537
BsrF I (2): 2094, 3369
BssH II (1): 4094
Bst1107 I (0)
BstB I (1): 654
BstE II (1): 3864
BstN I (13): 764, 932, 941, 2907, 2920, 3041, 3407, 3722, 4262, 4319, 4458, 4595, 4722
BstU I (26): 360, 1007, 1009, 1112, 1287, 1619, 2112, 2442, 3023, 3282, 3396, 3433, 3483, 3554, 3572, 3588, 3677, 3684, 3832, 3878, 3913, 4095, 4119, 4350, 4352, 4956
BstX I (3): 3478, 3607, 3730
BstY I (8): 914, 945, 1532, 1549, 2317, 2329, 2415, 2426
Bsu36 I (1): 4775
Cfr10 I (2): 2094, 3369
Cla I (0)
Dde I (8): 624, 1252, 1678, 2218, 2384, 2793, 4255, 4776
Dpn I (23): 660, 693, 915, 946, 998, 1497, 1533, 1550, 1808, 1854, 1872, 2213, 2318, 2330, 2408, 2416, 2427, 2502, 3594, 3698, 4071, 4673, 4769
Dpn II (23): 660, 693, 915, 946, 998, 1497, 1533, 1550, 1808, 1854, 1872, 2213, 2318, 2330, 2408, 2416, 2427, 2502, 3594, 3698, 4071, 4673, 4769
Dra I (5): 683, 792, 1599, 2291, 2310
Dra II (4): 289, 935, 936, 1198
Dra III (0)
Drd I (2): 1096, 2959
Dsa I (1): 4884
Dsa V (23): 764, 932, 941, 950, 959, 960, 1052, 1087, 1642, 1993, 2689, 2907, 2920, 3041, 3302, 3407, 3722, 4111, 4262, 4319, 4458, 4595, 4722
Eae I (5): 463, 975, 1786, 4357, 4560
Eag I (1): 975
Eam1105 I (1): 2174
Ear I (0)
Ear I (5): 338, 1381, 3185, 3302, 4654
Ecl136 II (0)
Eco47 III (0)
Eco57 I (1): 1492
Eco57 I (1): 2540
EcoN I (1): 264
EcoO109 I (4): 289, 935, 936, 1198
EcoR I (1): 954
EcoR II (13): 764, 932, 941, 2907, 2920, 3041, 3407, 3722, 4262, 4319, 4458, 4595, 4722
EcoR V (1): 4131

(continued on back)

Sheet 2 of 3 -Front

27-4597-01

Esp I (0)
Esp3 I (1): 4292
Esp3 I (1): 1051
Fnu4H I (30): 42, 974, 977, 1046, 1143, 1434, 1663, 1758, 1785, 1935, 2124, 2452, 2658, 2661, 2726, 2869, 3024, 3142, 3145, 3163, 3450, 3573, 3637, 3657, 3744, 4115, 4244, 4373, 4622, 4695
Fok I (7): 767, 1082, 1727, 2014, 2195, 3716, 3725
Fok I (3): 692, 887, 4615
Fsp I (2): 1954, 4692
Hae II (6): 2823, 3193, 3298, 4079, 4322, 4713
Hae III (22): 291, 464, 859, 937, 976, 1200, 1787, 2054, 2134, 2592, 3026, 3044, 3055, 3410, 3545, 3891, 4260, 4358, 4561, 4663, 4780, 4949
Hga I (3): 726, 1638, 4954
Hga I (8): 25, 1103, 2378, 2956, 3630, 3675, 3904, 3910
Hgi A I (5): 18, 1507, 1592, 2753, 3663
Hha I (34): 109, 138, 359, 1008, 1111, 1141, 1286, 1618, 1955, 2048, 2441, 2550, 2724, 2824, 2891, 3161, 3194, 3299, 3555, 3589, 3678, 3879, 3912, 4080, 4087, 4094, 4096, 4118, 4323, 4351, 4416, 4693, 4714, 4829
Hinc II (4): 182, 964, 1635, 4187
Hind II (4): 182, 964, 1635, 4187
Hind III (0)
Hinf I (10): 967, 2180, 2697, 3093, 3168, 3311, 4107, 4362, 4553, 4891
Hinp I (34): 109, 138, 359, 1008, 1111, 1141, 1286, 1618, 1955, 2048, 2441, 2550, 2724, 2824, 2891, 3161, 3194, 3299, 3555, 3589, 3678, 3879, 3912, 4080, 4087, 4094, 4096, 4118, 4323, 4351, 4416, 4693, 4714, 4829
Hpa I (0)
Hpa II (21): 951, 960, 1053, 1087, 1642, 1884, 1994, 2061, 2095, 2499, 2689, 2715, 2862, 3303, 3370, 3762, 4003, 4111, 4480, 4730, 4748
Hph I (12): 369, 438, 690, 1017, 1026, 1475, 2097, 2324, 3326, 3401, 4268, 4301
Hph I (3): 1453, 1694, 3859
Kas I (1): 4322
Kpn I (0)
Ksp632 I (0)
Ksp632 I (5): 338, 1381, 3185, 3302, 4654
Mae I (4): 270, 1986, 2321, 2574
Mae II (13): 1, 473, 800, 882, 1160, 1257, 1577, 1950, 2366, 3345, 3418, 4576, 4845
Mae III (19): 699, 985, 1062, 1157, 1519, 1707, 1860, 1918, 2249, 2532, 2648, 2711, 3342, 3865, 4581, 4601, 4821, 4847, 4905
Mam I (0)
Mbo I (23): 660, 693, 915, 946, 998, 1497, 1533, 1550, 1808, 1854, 1872, 2213, 2318, 2330, 2408, 2416, 2427, 2502, 3594, 3698, 4071, 4673, 4769
Mbo II (4): 412, 2411, 3664, 4764
Mbo II (13): 330, 342, 657, 1189, 1385, 1494, 1572, 2327, 3189, 3306, 3822, 4161, 4658
Mlu I (1): 3683
Mnl I (7): 1806, 2012, 2890, 3173, 4028, 4654, 4771
Mnl I (11): 903, 1004, 1034, 1202, 2149, 2230, 2630, 2954, 3535, 4342, 4800
Msc I (1): 463
Mse I (22): 190, 286, 445, 450, 684, 793, 1228, 1600, 1965, 2004, 2239, 2292, 2306, 2311, 2363, 3580, 3703, 3897, 4368, 4427, 4611, 4920
Msp I (21): 951, 960, 1053, 1087, 1642, 1884, 1994, 2061, 2095, 2499, 2689, 2715, 2862, 3303, 3370, 3762, 4003, 4111, 4480, 4730, 4748
Mun I (0)
Nae I (0)
Nar I (1): 4322
Nci I (10): 950, 959, 960, 1052, 1087, 1642, 1993, 2689, 3302, 4111
Nco I (0)
Nde I (0)
NgoM I (0)
Nhe I (0)
Nla III (22): 257, 497, 696, 716, 749, 856, 1043, 1148, 1235, 1340, 1733, 1769, 1847, 1857, 2348, 3068, 3292, 3819, 3998, 4019, 4171, 4543
Nla IV (18): 752, 935, 936, 937, 945, 1290, 1880, 2091, 2132, 2226, 2998, 3037, 3603, 3890, 4322, 4452, 4732, 4744
Not I (1): 974
Nru I (0)
Nsi I (0)
Nsp I (3): 496, 1042, 3067
NspB II (8): 1117, 1541, 2482, 2727, 3711, 4281, 4374, 4642
Pac I (0)
PaeR7 I (1): 969
PflM I (1): 3259
Ple I (2): 962, 2692
Ple I (4): 2180, 3168, 3311, 4107
Pme I (0)
Pml I (0)
PpuM I (0)
Psp1406 I (3): 1576, 1949, 3344
Pst I (1): 1933
Pvu I (2): 1807, 4672
Pvu II (3): 4281, 4374, 4642
Rsa I (3): 830, 1697, 3829
Rsr II (0)
Sac I (0)
Sac II (0)
Sal I (1): 964
Sau I (1): 4775
Sau3A I (23): 660, 693, 915, 946, 998, 1497, 1533, 1550, 1808, 1854, 1872, 2213, 2318, 2330, 2408, 2416, 2427, 2502, 3594, 3698, 4071, 4673, 4769
Sau96 I (15): 290, 752, 936, 937, 1199, 1815, 2037, 2054, 2133, 3545, 3890, 3891, 4235, 4259, 4663
Sca I (2): 829, 1696
Sce I (0)
Sce I (0)
ScrF I (23): 764, 932, 941, 950, 959, 960, 1052, 1087, 1642, 1993, 2689, 2907, 2920, 3041, 3302, 3407, 3722, 4111, 4262, 4319, 4458, 4595, 4722
Sdu I (7): 18, 936, 1507, 1592, 2753, 3663, 3890
Sec I (10): 932, 940, 941, 949, 959, 2907, 4318, 4457, 4594, 4884
SfaN I (5): 980, 1718, 3221, 3848, 4039
SfaN I (10): 759, 1074, 1479, 1919, 2971, 3191, 3717, 4046, 4804, 4817
Sfc I (4): 1933, 2611, 2802, 4940
Sfi I (0)
SgrA I (0)
Sin I (4): 752, 1815, 2037, 4235
Sma I (1): 959
SnaB I (0)
Spe I (0)
Sph I (0)
Spl I (0)
Spo I (0)
Srf I (0)
Sse8387 I (0)
Ssp I (2): 164, 1372
Stu I (0)
Sty I (0)
Swa I (1): 682
Taq I (11): 7, 308, 576, 655, 965, 970, 1525, 2969, 3262, 3621, 3643
Tfi I (4): 3093, 4362, 4553, 4891
Tth111 I (1): 1151
Xba I (0)
Xcm I (3): 3532, 4048, 4066
Xho I (1): 969
Xho II (8): 914, 945, 1532, 1549, 2317, 2329, 2415, 2426
Xma I (1): 959
Xmn I (2): 647, 1575

LOCUS pGEX-6P-1 4984 bp
BASE COUNT 1228 A 1207 C 1297 G 1252 T
ORIGIN

```
1  ACGTTATCGA CTGCACGGTG CACCAATGCT TCTGGCGTCA GGCAGCCATC GGAAGCTGTG
61  GTATGGCTGT GCAGGTCGTA AATCACTGCA TAATTCGTGT CGCTCAAGGC GCACTCCCGT
121 TCTGGATAAT GTTTTTTGCG CCGACATCAT AACGGTCTG GCAAATATTC TGAAATGAGC
181 TGTTGACAAT TAATCATCGG CTCGTATAAT GTGTGGAATT GTGAGCGGAT AACAATTTCA
241 CACAGGAAAC AGTATTCATG TCCCCTATAC TAGGTTATTG GAAAATTAAG GGCCTTGTGC
301 AACCCACTCG ACTTCTTTTG GAATATCTTG AAGAAAAATA TGAAGAGCAT TTGTATGAGC
361 GCGATGAAGG TGATAAATGG CGAAACAAAA AGTTTGAATT GGGTTTGGAG TTTCCCAATC
421 TTCCTTATTA TATTGATGGT GATGTTAAAT TAACACAGTC TATGGCCATC ATACGTTATA
481 TAGCTGACAA GCACAACATG TTGGGTGGTT GTCCAAAAGA GCGTGCAGAG ATTTCAATGC
541 TTGAAGGAGC GGTTTTGGAT ATTAGATACG GTGTTTCGAG AATTGCATAT AGTAAAGACT
601 TTGAAACTCT CAAAGTTGAT TTTCTTAGCA AGCTACCTGA AATGCTGAAA ATGTTCGAAG
661 ATCGTTTATG TCATAAAACA TATTTAAATG GTGATCATGT AACCCATCCT GACTTCATGT
721 TGTATGACGC TCTTGATGTT GTTTTATACA TGGACCCAAT GTGCCTGGAT GCGTTCCCAA
781 AATTAGTTTG TTTTAAAAAA CGTATTGAAG CTATCCCACA AATTGATAAG TACTTGAAAT
841 CCAGCAAGTA TATAGCATGG CCTTTCGAGG GCTGGCAAGC CACGTTTGGT GGTGGCGACC
901 ATCCTCCAAA ATCGGATCTG GAAGTCTGT TCCAGGGGCC CCTGGGATCC CCGGAATTCC
961 CGGGTCGACT CGAGCGGCCG CATCGTGACT GACTGACGAT CTGCCTCGCG CGTTTCGGTG
1021 ATGACGGTGA AAACCTCTGA CACATGCAGC TCCCGGAGAC GGTACACAGT TGTCTGTAAG
1081 CGGATGCCGG GAGCAGACAA GCCCGTCAGG GCGCGTCAGC GGGTGTGGC GGGTGTGGG
1141 GCGCAGCCAT GACCCAGTCA CGTAGCGATA GCGGAGTGTA TAATTCTTGA AGACGAAAGG
1201 GCCTCGTGAT ACGCCTATTT TTATAGGTTA ATGTCATGAT AATAATGGTT TCTTAGACGT
1261 CAGGTGGCAC TTTTCGGGGA AATGTGCGCG GAACCCCTAT TTGTTTATTT TTCTAAATAC
1321 ATTCAAATAT GTATCCGCTC ATGAGACAAT AACCTGATA AATGCTTCAA TAATATTGAA
1381 AAAGGAAGAG TATGAGTATT CAACATTTCC GTGTGCCCCT TATTCCTTTT TTTGCGGCAT
1441 TTTGCCTTCC TGTTTTTGCT CACCCAGAAA CGCTGGTGAA AGTAAAAGAT GCTGAAGATC
1501 AGTTGGGTGC ACGAGTGGGT TACATCGAAC TGGATCTCAA CAGCGGTAAG ATCCTTGAGA
1561 GTTTTCGCCC CGAAGAACGT TTTCCAATGA TGAGCACTTT TAAAGTTCTG CTATGTGGCG
1621 CGGTATTATC CCGTGTGAC GCCGGGCAAG AGCAACTCGG TCGCCGCATA CACTATTCTC
1681 AGAATGACTT GGTTGAGTAC TCACCAGTCA CAGAAAAGCA TCTTACGGAT GGCATGACAG
1741 TAAGAGAATT ATGCAGTGCT GCCATAACCA TGAGTGATAA CACTGCGGCC AACTTACTTC
1801 TGACAACGAT CGGAGGACCG AAGGAGCTAA CCGCTTTTTT GCACAACATG GGGGATCATG
1861 TAACTCGCCT TGATCGTTGG GAACCGGAGC TGAATGAAGC CATAACAAAC GACGAGCGTG
1921 ACACCACGAT GCCTGCAGCA ATGGCAACAA CGTTGCGCAA ACTATTAACT GGCGAACTAC
1981 TTA CTAGC TTCCC GGCAA CAATTAATAG ACTGGATGGA GCGGATAAAA GTTGCAGGAC
2041 CACTTCTGCG CTCGGCCCTT CCGGCTGGCT GGTTTATTGC TGATAAATCT GGAGCCGGTG
2101 AGCGTGGGTC TCGCGGTATC ATTGCAGCAC TGGGGCCAGA TGGTAAGCCC TCCCGTATCG
2161 TAGTTATCTA CACGACGGGG AGTCAGGCAA CTATGGATGA ACGAAATAGA CAGATCGCTG
2221 AGATAGGTGC CTCACTGATT AAGCATTGGT AACTGTCAGA CCAAGTTTAC TCATATATAC
2281 TTTAGATTGA TTTAAAACCTT CATTTTTAAT TTAAAAGGAT CTAGGTGAAG ATCCTTTTTG
2341 ATAATCTCAT GACCAAAATC CCTTAACGTG AGTTTTCGTT CCACTGAGCG TCAGACCCCG
2401 TAGAAAAGAT CAAAGGATCT TCTTGAGATC CTTTTTTTCT GCGCGTAATC TGCTGCTTGC
2461 AAACAAAAAA ACCACCGCTA CCAGCGGTGG TTTGTTTGCC GGATCAAGAG CTACCAACTC
2521 TTTTTCGGAA GGTAACTGGC TTCAGCAGAG CGCAGATACC AAATACTGTC CTTCTAGTGT
2581 AGCCGTAGTT AGGCCACCAC TTCAAGAACT CTGTAGCACC GCCTACATAC CTCGCTCTGC
2641 TAATCCTGTT ACCAGTGGCT GCTGCCAGTG GCGATAAGTC GTGTCTTACC GGGTTGGACT
2701 CAAGACGATA GTTACCGGAT AAGGCGCAGC GGTGCGGCTG AACGGGGGGT TCGTGCACAC
2761 AGCCCAGCTT GGAGCGAACG ACCTACACCG AACTGAGATA CCTACAGCGT GAGCTATGAG
2821 AAAGCGCCAC GCTTCCCGAA GGGAGAAAGG CGGACAGGTA TCCGGTAAGC GGCAGGGTCG
2881 GAACAGGAGA GCGCACGAGG GAGCTTCCAG GGGGAAACGC CTGGTATCTT TATAGTCTTG
2941 TCGGGTTTCG CCACCTCTGA CTTGAGCGTC GATTTTTGTG ATGCTCGTCA GGGGGGCGGA
3001 GCCTATGGAA AAACGCCAGC AACGCGGCCT TTTTACGGTT CCTGGCCTTT TGCTGGCCTT
3061 TTGCTCACAT GTTCTTTCTT CCGTTATCCC CTGATTCTGT GGATAACCGT ATTACCGCCT
3121 TTGAGTGAGC TGATACCGCT CGCCGACGCC GAACGACCGA GCGCAGCGAG TCAGTGAGCG
```

3181 AGGAAGCGGA AGAGCGCCTG ATGCGGTATT TTCTCCTTAC GCATCTGTGC GGTATTTAC
3241 ACCGCATAAA TTCCGACACC ATCGAATGGT GCAAAACCTT TCGCGGTATG GCATGATAGC
3301 GCCCGBAAGA GAGTCAATTC AGGGTGGTGA ATGTGAAACC AGTAACGTTA TACGATGTGC
3361 CAGAGTATGC CGGTGTCTCT TATCAGACCG TTTCCCGCGT GGTGAACCAG GCCAGCCACG
3421 TTTCTGCGAA AACCGGGGAA AAAGTGAAG CGGCGATGGC GGAGCTGAAT TACATTCCCA
3481 ACCGCGTGGC ACAACAACCTG GCGGGCAAAC AGTCGTTGCT GATTGGCGTT GCCACCTCCA
3541 GTCTGGCCCT GCACGCGCCG TCGCAAATTG TCGCGGCGAT TAAATCTCGC GCCGATCAAC
3601 TGGGTGCCAG CGTGGTGGTG TCGATGGTAG AACGAAGCGG CGTCGAAGCC TGTAAGCGG
3661 CGGTGCACAA TCTTCTCGCG CAACGCGTCA GTGGGCTGAT CATTAACTAT CCGCTGGATG
3721 ACCAGGATGC CATTGCTGTG GAAGCTGCCT GCACTAATGT TCCGGCGTTA TTTCTTGATG
3781 TCTCTGACCA GACACCCATC AACAGTATTA TTTTCTCCCA TGAAGACGGT ACGCGACTGG
3841 GCGTGGAGCA TCTGGTCGCA TTGGGTACC AGCAAATCGC GCTGTTAGCG GGCCCATTA
3901 GTTCTGTCTC GCGCGTCTG CGTCTGGCTG GCTGGCATAA ATATCTCACT CGCAATCAAA
3961 TTCAGCCGAT AGCGGAACGG GAAGGCGACT GGAGTGCCAT GTCCGGTTTT CAACAAACCA
4021 TGCAAATGCT GAATGAGGGC ATCGTTCCCA CTGCGATGCT GGTGCCAAC GATCAGATGG
4081 CGCTGGGCGC AATGCGCGCC ATTACCGAGT CCGGGCTGCG CGTTGGTGCG GATATCTCGG
4141 TAGTGGGATA CGACGATACC GAAGACAGCT CATGTTATAT CCCGCCGTCA ACCACCATCA
4201 AACAGGATTT TCGCCTGCTG GGGCAAACCA GCGTGGACCG CTTGCTGCAA CTCTCTCAGG
4261 GCCAGGCGGT GAAGGGCAAT CAGCTGTTGC CCGTCTCACT GGTGAAAAGA AAAACCACCC
4321 TGGCGCCCAA TACGCAAACC GCCTCTCCCC GCGCGTTGGC CGATTCATTA ATGCAGCTGG
4381 CACGACAGGT TTCCCGACTG GAAAGCGGGC AGTGAGCGCA ACGCAATTAA TGTGAGTTAG
4441 CTCACTCATT AGGCACCCCA GGCTTTACAC TTTATGCTTC CGGCTCGTAT GTTGTGTGGA
4501 ATTGTGAGCG GATAACAATT TCACACAGGA AACAGCTATG ACCATGATTA CGGATTCACT
4561 GGCCGTCGTT TTACAACGTC GTGACTGGGA AAACCCTGGC GTTACCCAAC TTAATCGCCT
4621 TGCAGCACAT CCCCCTTTTCG CCAGCTGGCG TAATAGCGAA GAGGCCCGCA CCGATCGCCC
4681 TTCCCAACAG TTGCGCAGCC TGAATGGCGA ATGGCGCTTT GCCTGGTTTC CGGCACCAGA
4741 AGCGGTGCCG GAAAGCTGGC TGGAGTGCGA TCTTCCTGAG GCCGATACTG TCGTCGTCCC
4801 CTCAAACTGG CAGATGCACG GTTACGATGC GCCCATCTAC ACCAACGTAA CCTATCCCAT
4861 TACGGTCAAT CCGCCGTTTG TTCCACGGA GAATCCGACG GGTGTTACT CGCTCACATT
4921 TAATGTTGAT GAAAGCTGGC TACAGGAAGG CCAGACGCGA ATTATTTTTG ATGGCGTTGG
4981 AATT