

	IPG Strips					IPG Buffers							
	Strip length					pH range							
	24 cm	18 cm	13 cm	11 cm	7 cm	3.5-5.0	5.5-6.7	4-7	6-11	7-11 NL	3-10 NL	3-10	3-11 NL
Narrow													
3.5-4.5	×		×	×	×	●							
5.3-6.5	×	×	×	×	×		●						
6.2-7.5	×	×	×	×	×				●				
Medium													
3-5.6 NL	×	×	×	×	×	●							
3-7 NL	×							●					
4-7	×	×	×	×	×			●					
6-9	×	×							●				
6-11		×	×	×	×				●				
7-11 NL	×	×	×	×	×					●			
Wide													
3-10	×	×	×	×	×							●	
3-11 NL	×	×	×	×	×								●
3-10 NL	×	×	×		×						●		